



IBM iConnect[®] Enterprise Archive 13.1

THE INTEGRATION STATEMENT

Licensed Materials - Property of IBM

© Copyright IBM Corporation. 1998, 2020.

IBM, the IBM logo, IBM Watson, and Watson Health are trademarks of IBM Corporation in the United States, other countries or both. The use of IBM trademarks herein are by authorization and license from the IBM Corporation.

Microsoft, Windows, Windows NT, and the Windows logo are trademarks of Microsoft Corporation in the United States, other countries, or both.

All other company or product names are trademarks of their respective companies.

Licensees of this program who wish to have information about it for the purpose of enabling: (i) the exchange of information between independently created programs and other programs (including this one) and (ii) the mutual use of the information which has been exchanged, should contact the IBM Director of Licensing. Such information may be available, subject to appropriate terms and conditions, including in some cases, payment of a fee.

IBM Director of Licensing
IBM Corporation
North Castle Drive, MD-NC119
Armonk, NY 10504-1785

IBM grants limited permission to licensees to make hard copy or other reproductions of any machine-readable documentation, provided that each such reproduction shall carry the IBM copyright notices and that use of the reproduction shall be governed by the terms and conditions specified by IBM in the Licensed Program Specifications. Any reproduction or use beyond the limited permission granted herein shall be a breach of the license agreement and an infringement of the applicable copyrights.

Materials displayed or reproduced by this program may be protected by copyright and/or contract restrictions of IBM and/or others. The user is responsible for having permission to display and reproduce such materials and for including applicable copyright notices and legends.

U.S. GOVERNMENT RESTRICTED RIGHTS:

Use, duplication, or disclosure restricted by GSA ADP Schedule Contract with IBM Corporation.

INSTRUCTIONS FOR USE:


Instructions for use are available electronically in PDF format. Paper format can be requested for supported languages by contacting Customer Support, and will be provided free of charge within seven calendar days of receiving the request. The paper format can also be provided at the time of delivery if so requested at the time of order.

INTENDED USE:

The IBM iConnect Enterprise Archive is a software system that supports storage and communications of medical images and data. Qualified system administrators monitor and maintain the software in their healthcare enterprise environment.

CAUTION: U.S. federal law restricts this device to sale by, or on the order of, a physician.

CANADIAN DEVICE IDENTIFIER: 100-00002-00

 Merge Healthcare Incorporated
900 Walnut Ridge Drive
Hartland, WI 53029

The symbols glossary is provided electronically at <https://developer.ibm.com/watsonhealth/ibm-merge-statements-and-patches/>.

For application support or to report issues with user documentation, contact Customer Support:

- In North America: call toll-free 1-877-741-5369
- Outside of North America: +31.40.299.0773
- Email: WHISupport@us.ibm.com

Part	Date	Version	Description
ECM-22997	March 2020	1.0	Initial version for 13.1

The latest revision of this document can be found in [IBM Watson Health Community \(https://mergecustomer.force.com/mergeusercommunity/login\)](https://mergecustomer.force.com/mergeusercommunity/login).



IHE Integration Statement

The following table outlines the IHE Integration Statement.

IHE Integration Statement		Date
		2020
Vendor	Product Name	Version
Merge Healthcare Incorporated	IBM iConnect® Enterprise Archive	13.1
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Patient Information Reconciliation	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Consistent Presentation of Images	Image Manager / Image Archive	None
Presentation of Grouped Procedures	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Evidence Documents	Image Manager / Image Archive	None
Key Image Note	Image Manager / Image Archive	None
Simple Image and Numeric Report	Report Repository	None
	Enterprise Report Repository	None
	External Report Repository Access	None

Access to Radiology Information	Image Manager / Image Archive	None
	Report Repository	None
	External Report Repository Access	None
Nuclear Medicine	Image Manager / Image Archive	None
Cardiac Catheterization Workflow	Image Manager / Image Archive	Intermittently Connected Modality, Cardiac Cath
	Performed Procedure Step Manager	None
Echo Workflow	Image Manager / Image Archive	Echocardiography, Intermittently Connected Modality
	Performed Procedure Step Manager	None
Evidence Documents - Cardiology	Image Manager / Image Archive	None
Consistent Time (CT)	Time Client	None
Mammography	Image Manager / Image Archive	None
Stress Testing Workflow	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Image Document Source	Set of DICOM Instances
	Document Repository	None
Cross-Enterprise Document Sharing (XDS.b)	Document Registry	Patient Identity Feed (HL7 V2), Patient Identity Feed (HL7 V3), On-Demand Documents, Reference ID, Document Metadata Update
	Document Repository	None
	Integrated Document Source/Repository	Document Replace
Audit Trail and Node Authentication (ATNA)	Secure Application	None
Patient Identifier Cross-reference (PIX)	Patient Identity Cross-reference Manager	None
	Patient Identifier Cross-reference Consumer	PIX Update Notification

Patient Demographics Query (PDQ)	Patient Demographics Supplier	None
	Patient Demographics Consumer	None
Imaging Object Change Management (IOCM)	Image Manager/Archive	None
	Change Requester	None
Document Metadata Subscription (DSUB)	Document Metadata Notification Broker	Patient-Independent Subscription, Folder Subscription
	Document Metadata Publisher	None
Multi-Patient Queries (MPQ)	Document Registry	None
Web-based Image Capture (WIC)	Image Manager	None
Internet address for vendor's IHE information:	www.merge.com/Support/Resources/IHE-Statements.aspx	
Links to Standards Conformance Statements for the Implementation		
HL7	Refer to Merge Communities -- http://www.merge.com/Support/Support-Portal.aspx	
DICOM	www.merge.com/Support/Resources/DICOM-Conformance-Statements.aspx	
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org